



The **BEST** for  
Life Science

## TransStart® TopTaq DNA Polymerase

Cat. No. AP151

Concentration 2.5 units/μl

Storage at -20°C for two years

### Description

TransStart® TopTaq DNA Polymerase is an engineered version of *Taq* DNA Polymerase combined with TransStart® technique. One binding protein binds to double-strand DNA template, preventing polymerase activity at room temperature. Other two binding proteins bind primers, preventing primer-dimer formation. Blocking proteins are released from primers and templates during the initial denaturation. This double blocking method has higher efficiency than antibody based, or chemically modified hot start PCR.

### Highlights

- Compared with TransStart® *Taq* DNA Polymerase, TransStart® TopTaq DNA Polymerase has higher amplification efficiency, specificity and sensitivity.
- TransStart® TopTaq DNA Polymerase offers 18-fold fidelity as compared to EasyTaq® DNA Polymerase.
- The specificity is higher than antibody based or chemically modified hot start DNA polymerases.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be directly cloned into pEASY®-T vectors.
- Reduced nonspecific amplification and primer dimer formation.
- Different from *Taq* antibody, no risk of contamination from mammalian DNA.
- Different from chemical modification, long denaturing step is not needed.
- Amplification of genomic DNA fragment up to 15 kb.

### Applications

- Complex templates
- GC/AT-rich templates
- Multiplex PCR
- High yield PCR

### Unit Definition

One unit of TransStart® TopTaq DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

TransStart® TopTaq DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of TransStart® TopTaq DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

### Storage Buffer

20 mM Tris-HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, 50% glycerol, stabilizers

### 10×TransStart® TopTaq Buffer with 20 mM MgSO<sub>4</sub>

500 mM Tris-HCl (pH 9.0), 200 mM (NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>, 20 mM MgSO<sub>4</sub>, others

### GC Enhancer

For better amplification of GC rich or complex templates, we recommend adding GC enhancer to PCR reaction. GC enhancer is provided at 10× concentration and can be used at 0.5×-5× concentration.

### Kit Contents

Component	AP151-01/11	AP151-02/12	AP151-03/13
<i>TransStart</i> <sup>®</sup> <i>TopTaq</i> DNA Polymerase	250 U×1	500 U×1	500 U×6
10× <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> Buffer	1.2 ml	1.2 ml×2	1.2 ml×12
2.5 mM dNTPs	-/800 µl×1	-/800 µl×2	-/1.2 ml×8
10×GC Enhancer	200 µl×1	400 µl×1	1 ml×1
6×DNA Loading Buffer	500 µl×1	1 ml×1	1 ml×2

### Reaction Components

Component	Volume	Final Concentration
Template DNA	Variable	as required
Forward Primer (10 µM)	1 µl	0.2 µM
Reverse Primer (10 µM)	1 µl	0.2 µM
10× <i>TransStart</i> <sup>®</sup> <i>TopTaq</i> Buffer	5 µl	1×
2.5 mM dNTPs	4 µl	0.2 mM
<i>TransStart</i> <sup>®</sup> <i>TopTaq</i> DNA Polymerase	0.5-1 µl	1.25-2.5 units
Nuclease-free Water	Variable	-
Total volume	50 µl	-

### Thermal cycling conditions

94°C	2-5 min	} 30-35 cycles
94°C	30 sec	
50-60°C	30 sec	
72°C	1-2 kb/min	
72°C	5-10 min	

### Note

- A final concentration of 2 mM MgSO<sub>4</sub> is sufficient for most targets amplification. For some targets, more Mg<sup>2+</sup> may be required.
- For optimal results, we recommend to use the 100 mM MgSO<sub>4</sub> stock to prepare a titration from 2 mM to 4 mM (final concentration) in 0.25 mM increments.
- 0.5 µl (2.5 units) enzyme is enough for per 50 µl reaction. For better amplification, up to 1 µl (5 units) enzyme can be used.

FOR RESEARCH USE ONLY