

Pre-stained protein Marker (10~180 kDa)

Cat. No: E-BC-R273

Size: 250 µL/ 500 µL/ 2.5 mL

Cat	Products	250 µL	500 µL	2.5 mL	Storage
E-BC-R273	Pre-stained Protein Marker (10~180 kDa)	250 µL	250 µL× 2	250 µL× 10	-20 °C

Introduction

This product consists of 10 pre-stained proteins with a molecular weight range of 10~180 kDa, 10 kDa is a green band, 72 kDa is an orange band, and the rest are blue bands. It can be used to directly observe the protein electrophoresis status and clearly judge the transmembrane effect of Western Blot.

Instructions

1. Thaw the product at room temperature. Do not heat it.
2. Make sure that the solution is fully mixed before use.
3. Add an appropriate volume of pre-stained Marker to the well of SDS-PAGE gel. It is recommended to add 5 µL of Pre-stained protein Marker to the microgel of 1 mm thick mini-gel and 10 µL of Pre-stained protein Marker to the full length gel.
4. Please store at -20 °C in time after use.

Storage Buffer


62.5 mM Tris•H₃PO₄ (pH7.5 at 25 °C), 1 mM EDTA, 2% (w/v) SDS, 10 mM DTT, 1 mM NaN₃ and 33% (v/v) Glycerin.

Storage

This product can be stored at -20 °C for 12 months or at 2~8 °C for 3 months.

Caution

Do not boil this product.



For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about products, please visit www.elabscience.com

Focus on lab research
 Service for life science