

MEM Medium Ingredient Formula_Table 1

Ingredient(mg/L) / Cat.	PM150410	PM150410A	PM150412	PM150412A	PM150414	PM150415	PM150416	PM150417
Inorganic Salt								
CaCl2 (Anhydrous)	200	200	200	200	200	200	200	200
MgSO4 (Anhydrous)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO3	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6800	6800	6800	6800	6350	6350	6800	6350
NaH2PO4 (Anhydrous)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	–	–	–	–	–	–	434	434
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	292	292	–	–	292	–	–	–
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
L-Alanine	8.9	8.9	8.9	8.9	8.9	8.9	8.9	8.9
L -Asparagine	15	15	15	15	15	15	15	15
L-Aspartic Acid	13.3	13.3	13.3	13.3	13.3	13.3	13.3	13.3
L-Glutamic Acid	14.7	14.7	14.7	14.7	14.7	14.7	14.7	14.7
L-Proline	11.5	11.5	11.5	11.5	11.5	11.5	11.5	11.5
L-Serine	10.5	10.5	10.5	10.5	10.5	10.5	10.5	10.5
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	1000	1000	1000	1000	1000	1000	1000	1000
HEPES	–	–	–	–	5958	5958	–	5958
Phenol Red	10	10	10	10	10	10	10	10
Penicillin G Sodium Salt	–	100kU	–	100kU	–	–	–	–
Streptomycin Sulfate	–	100	–	100	–	–	–	–

MEM Medium Ingredient Formula_Table 2

Ingredient(mg/L) / Cat.	PM150418	PM150419	PM150420	PM150430	PM150431	PM150432	PM150411	PM150411A
Inorganic Salt								
CaCl2 (Anhydrous)	200	200	200	200	200	200	200	200
MgSO4 (Anhydrous)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO3	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6800	6800	6350	6350	6800	6350	6800	6800
NaH2PO4 (Anhydrous)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	-	-	-	-	434	434	-	-
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	292	-	292	-	-	-	292	292
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	7.5	7.5	7.5	7.5	7.5	7.5	-	-
L-Alanine	8.9	8.9	8.9	8.9	8.9	8.9	-	-
L-Asparagine	15	15	15	15	15	15	-	-
L-Aspartic Acid	13.3	13.3	13.3	13.3	13.3	13.3	-	-
L-Glutamic Acid	14.7	14.7	14.7	14.7	14.7	14.7	-	-
L-Proline	11.5	11.5	11.5	11.5	11.5	11.5	-	-
L-Serine	10.5	10.5	10.5	10.5	10.5	10.5	-	-
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	1000	1000	1000	1000	1000	1000	1000	1000
HEPES	-	-	5958	5958	-	5958	-	-
Phenol Red	-	-	-	-	-	-	10	10
Penicillin G Sodium Salt	-	-	-	-	-	-	-	100kU
Streptomycin Sulfate	-	-	-	-	-	-	-	100

MEM Medium Ingredient Formula_Table 3

Ingredient(mg/L) / Cat.	PM150413	PM150413A	PM150433	PM150434	PM150435	PM150436	PM150437	PM150438
Inorganic Salt								
CaCl ₂ (无水)	200	200	200	200	200	200	200	200
MgSO ₄ (无水)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO ₃	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6800	6800	6350	6350	6800	6350	6800	6800
NaH ₂ PO ₄ (无水)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	-	-	-	-	434	434	-	-
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	-	-	292	-	-	-	292	-
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	-	-	-	-	-	-	-	-
L-Alanine	-	-	-	-	-	-	-	-
L-Asparagine	-	-	-	-	-	-	-	-
L-Aspartic Acid	-	-	-	-	-	-	-	-
L-Glutamic Acid	-	-	-	-	-	-	-	-
L-Proline	-	-	-	-	-	-	-	-
L-Serine	-	-	-	-	-	-	-	-
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	1000	1000	1000	1000	1000	1000	1000	1000
HEPES	-	-	5958	5958	-	5958	-	-
Phenol Red	10	10	10	10	10	10	-	-
Penicillin G Sodium Salt	-	100kU	-	-	-	-	-	-
Streptomycin Sulfate	-	100	-	-	-	-	-	-

MEM Medium Ingredient Formula_Table 4

Ingredient(mg/L) / Cat. No.	PM150439	PM150440	PM150441	PM150442	PM150443	PM150444	PM150445	PM150446
Inorganic Salt								
CaCl2 (Anhydrous)	200	200	200	200	200	200	200	200
MgSO4 (Anhydrous)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO3	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6350	6350	6800	6350	6800	6800	6350	6350
NaH2PO4 (Anhydrous)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	-	-	434	434	-	-	-	-
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	292	-	-	-	292	-	292	-
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	-	-	-	-	7.5	7.5	7.5	7.5
L-Alanine	-	-	-	-	8.9	8.9	8.9	8.9
L-Asparagine Monohydrate	-	-	-	-	15	15	15	15
L-Aspartic Acid	-	-	-	-	13.3	13.3	13.3	13.3
L-Glutamic Acid	-	-	-	-	14.7	14.7	14.7	14.7
L-Proline	-	-	-	-	11.5	11.5	11.5	11.5
L-Serine	-	-	-	-	10.5	10.5	10.5	10.5
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	1000	1000	1000	1000	-	-	-	-
HEPES	5958	5958	5958	5958	-	-	5958	5958
Phenol Red	-	-	-	-	10	10	10	10
Penicillin G Sodium Salt	-	-	-	-	-	-	-	-
Streptomycin Sulfate	-	-	-	-	-	-	-	-

MEM Medium Ingredient Formula_Table 5

Ingredient(mg/L) / Cat.	PM150447	PM150448	PM150449	PM150450	PM150451	PM150452	PM150453	PM150454
Inorganic Salt								
CaCl2 (Anhydrous)	200	200	200	200	200	200	200	200
MgSO4 (Anhydrous)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO3	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6800	6350	6800	6800	6350	6350	6800	6350
NaH2PO4 (Anhydrous)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	434	434	–	–	–	–	434	434
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	–	–	292	–	292	–	–	–
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
L-Alanine	8.9	8.9	8.9	8.9	8.9	8.9	8.9	8.9
L-Asparagine	15	15	15	15	15	15	15	15
L-Aspartic Acid	13.3	13.3	13.3	13.3	13.3	13.3	13.3	13.3
L-Glutamic Acid	14.7	14.7	14.7	14.7	14.7	14.7	14.7	14.7
L-Proline	11.5	11.5	11.5	11.5	11.5	11.5	11.5	11.5
L-Serine	10.5	10.5	10.5	10.5	10.5	10.5	10.5	10.5
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	–	–	–	–	–	–	–	–
HEPES	–	5958	–	–	5958	5958	–	5958
Phenol Red	10	10	–	–	–	–	–	–
Penicillin G Sodium Salt	–	–	–	–	–	–	–	–
Streptomycin Sulfate	–	–	–	–	–	–	–	–

MEM Medium Ingredient Formula_Table 6

Ingredient(mg/L) / Cat.	PM150455	PM150456	PM150457	PM150458	PM150459	PM150460	PM150461	PM150462
Inorganic Salt								
CaCl2 (Anhydrous)	200	200	200	200	200	200	200	200
MgSO4 (Anhydrous)	97.67	97.67	97.67	97.67	97.67	97.67	97.67	97.67
KCl	400	400	400	400	400	400	400	400
NaHCO3	2200	2200	2200	2200	2200	2200	2200	2200
NaCl	6800	6800	6350	6350	6800	6350	6800	6800
NaH2PO4 (Anhydrous)	122	122	122	122	122	122	122	122
Amino Acid								
L-Alanyl-L-Glutamine	-	-	-	-	434	434	-	-
L-Arginine Hydrochloride	126	126	126	126	126	126	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	31.3	31.3	31.3
L-Glutamine	292	-	292	-	-	-	292	-
L-Histidine Hydrochloride	42	42	42	42	42	42	42	42
L-Isoleucine	52	52	52	52	52	52	52	52
L-Leucine	52	52	52	52	52	52	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	72.5	72.5	72.5
L-Methionine	15	15	15	15	15	15	15	15
L-Phenylalanine	32	32	32	32	32	32	32	32
L-Threonine	48	48	48	48	48	48	48	48
L-Tryptophan	10	10	10	10	10	10	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	51.9	51.9	51.9
L-Valine	46	46	46	46	46	46	46	46
Glycine	-	-	-	-	-	-	-	-
L-Alanine	-	-	-	-	-	-	-	-
L-Asparagine	-	-	-	-	-	-	-	-
L-Aspartic Acid	-	-	-	-	-	-	-	-
L-Glutamic Acid	-	-	-	-	-	-	-	-
L-Proline	-	-	-	-	-	-	-	-
L-Serine	-	-	-	-	-	-	-	-
Vitamin								
Choline Chloride	1	1	1	1	1	1	1	1
D-Calcium Pantothenate	1	1	1	1	1	1	1	1
Folic Acid	1	1	1	1	1	1	1	1
Niacinamide	1	1	1	1	1	1	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	1	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	1	1	1
Meso-inositol	2	2	2	2	2	2	2	2
Others								
D-Glucose	-	-	-	-	-	-	-	-
HEPES	-	-	5958	5958	-	5958	-	-
Phenol Red	10	10	10	10	10	10	-	-
Penicillin G Sodium Salt	-	-	-	-	-	-	-	-
Streptomycin Sulfate	-	-	-	-	-	-	-	-

MEM Medium Ingredient Formula_Table 7

Ingredient(mg/L) / Cat.	PM150463	PM150464	PM150465	PM150466	PM150467	PM150468	PM150469	PM150470
Inorganic Salt								
CaCl ₂ (Anhydrous)	200	200	200	200	200	2000	200	185
MgSO ₄ (Anhydrous)	97.67	97.67	97.67	97.67	97.67	976.7	97.67	97.67
KCl	400	400	400	400	400	4000	400	400
KH ₂ PO ₄	-	-	-	-	-	-	-	60
NaHCO ₃	2200	2200	2200	2200	1500	-	3700	350
NaCl	6350	6350	6800	6350	6800	68000	6800	8000
NaH ₂ PO ₄ (Anhydrous)	122	122	122	122	122	1220	122	122
Na ₂ HPO ₄ (Anhydrous)	-	-	-	-	-	-	-	47.88
Amino Acid								
L-Alanyl-L-Glutamine	-	-	434	434	-	-	-	-
L-Arginine Hydrochloride	126	126	126	126	126	1260	126	126
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	31.3	3128	31.3	31.3
L-Glutamine	292	-	-	-	292	-	-	292
L-Histidine Hydrochloride	42	42	42	42	42	420	42	42
L-Isoleucine	52	52	52	52	52	520	52	52
L-Leucine	52	52	52	52	52	520	52	52
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	72.5	725	72.5	72.5
L-Methionine	15	15	15	15	15	150	15	15
L-Phenylalanine	32	32	32	32	32	320	32	32
L-Threonine	48	48	48	48	48	480	48	48
L-Tryptophan	10	10	10	10	10	100	10	10
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	51.9	519	51.9	51.9
L-Valine	46	46	46	46	46	460	46	46
Glycine	-	-	-	-	7.5	-	-	-
L-Alanine	-	-	-	-	8.9	-	-	-
L-Asparagine	-	-	-	-	15	-	-	-
L-Aspartic Acid	-	-	-	-	13.3	-	-	-
L-Glutamic Acid	-	-	-	-	14.7	-	-	-
L-Proline	-	-	-	-	11.5	-	-	-
L-Serine	-	-	-	-	10.5	-	-	-
Vitamin								
Choline Chloride	1	1	1	1	1	10	1	1
D-Calcium Pantothenate	1	1	1	1	1	10	1	1
Folic Acid	1	1	1	1	1	10	1	1
Niacinamide	1	1	1	1	1	10	1	1
Pyridoxal Hydrochloride	1	1	1	1	1	10	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.1	1	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	1	10	1	1
Meso-inositol	2	2	2	2	2	20	2	2
Others								
D-Glucose	-	-	-	-	-	10000	6000	1000
HEPES	5958	5958	-	5958	-	-	-	-
Phenol Red	-	-	-	-	10	100	10	10
Sodium Pyruvate	-	-	-	-	110	-	-	-

MEM Medium Ingredient Formula_Table 8

Ingredient(mg/L) / Cat.	PM150471	PM150472	PM150473	PM150474	PM150475	PM150476	PM150477	PM150462
Inorganic Salt								
CaCl ₂ (Anhydrous)	–	200	200	200	285	200	–	200
FeCl ₃ •6H ₂ O	–	0.54	0.54	0.54	–	–	–	–
Fe(NO ₃) ₃ •9H ₂ O	–	–	–	–	1	0.1	–	–
MgCl ₂ (Anhydrous)	–	85.75	85.75	85.75	–	–	93.68	–
MgSO ₄ (Anhydrous)	97.67	12.22	12.22	12.22	195	97.67	–	97.67
KCl	400	400	400	400	800	400	400	400
NaHCO ₃	2200	2200	2200	2200	8400	2750	2000	850
NaCl	6800	6800	6800	6800	12800	6400	6500	6800
NaH ₂ PO ₄ (Anhydrous)	122	130	130	130	217.4	107.8	1154	122
ZnSO ₄ •6H ₂ O	–	0.14	0.14	0.14	–	–	–	–
Amino Acid								
L-Alanyl-L-Glutamine	–	–	–	–	–	–	–	–
L-Arginine Hydrochloride	126	126	126	126	504	42	126	126.4
L-Cystine Dihydrochloride	31.3	31.3	31.3	31.3	125.2	31.3	32.4	24.02
L-Glutamine	–	292	292	292	870	292	292	–
L-Histidine Hydrochloride	42	42	42	42	168	21	42	41.92
L-Isoleucine	52	52	52	52	208	52	52	52.4
L-Leucine	52	131	131	131	208	52	52	52.4
L-Lysine Hydrochloride	72.5	72.5	72.5	72.5	290	73.1	72.5	73.06
L-Methionine	15	15	15	15	60	15	15	14.92
L-Phenylalanine	32	32	32	32	128	33	32	33.02
L-Threonine	48	48	48	48	192	47.6	48	47.64
L-Tryptophan	10	10	10	10	40	8	10	10.2
L-Tyrosine Disodium Salt	51.9	51.9	51.9	51.9	207.6	52.19	54.52	36.22
L-Valine	46	46	46	46	155	46.8	46	46.86
Glycine	–	–	–	–	–	–	–	10.5
L-Alanine	–	–	–	–	–	–	–	8.9
L-Asparagine	–	60	60	60	–	–	–	15
L-Aspartic Acid	–	–	–	–	–	–	–	13.3
L-Glutamic Acid	–	–	–	–	–	–	–	14.7
L-Proline	–	–	–	–	–	–	–	11.5
L-Serine	–	42	42	42	–	–	–	10.5
Vitamin								
Biotin	–	0.1	0.1	0.1	–	–	–	–
Choline Chloride	1	56	56	56	4	2	1	1
D-Calcium Pantothenate	1	1	1	1	4	2	1	1
Folic Acid	1	2.2	2.2	2.2	4	2	1	1
Niacinamide	1	1	1	1	4	2	1	1
Pyridoxal Hydrochloride	1	1	1	1	4	2	1	1
Riboflavin	0.1	0.1	0.1	0.1	0.4	0.2	0.1	0.1
Thiamine Hydrochlorid	1	1	1	1	4	2	1	1
Vitamin B12	–	1.36	1.36	1.36	–	–	–	–
Meso-inositol	2	36	36	36	8	3.6	2	2
Others								
D-Glucose	1000	2000	2000	2000	9000	4500	2000	1000
Linoleic Acid	–	0.084	0.084	0.084	–	–	–	–
Znso ₄	–	0.2	0.2	0.2	–	–	–	–
Putrescine	–	0.16	0.16	0.16	–	–	–	–
Sodium Pyruvate	–	110	110	110	–	–	–	–
Phenol Red	10	10	–	10	–	15	10	10
Gentamicin Sulfate	–	–	0.05	–	–	–	–	–