

Hank's Solution (HBSS) Ingredient Formula

Ingredient (mg/L) / Cat. No.	PB180321	PB180322	PB180323	PB180324
Inorganic Salt				
CaCl ₂ (Anhydrous)	–	–	140	140
MgCl ₂ •6H ₂ O	–	–	100	100
MgSO ₄ •7H ₂ O	–	–	100	100
KCl	400	400	400	400
KH ₂ PO ₄	60	60	60	60
NaHCO ₃	350	350	350	350
NaCl	8000	8000	8000	8000
Na ₂ HPO ₄ •12H ₂ O	121	121	121	121
Others				
D-Glucose	1000	1000	1000	1000
Phenol Red	–	10	–	10