

Freezing Medium (Serum-free & animal origin-free)

Cat. No: PB180438
Size: 100mL

General Information

Product Form	Liquid
Concentration	1×
Storage	2-8°C. If not used for a long time, please store at -20°C.
Shipping	Ice bag
Expiration date	36 months

Background

Freezing Medium (Serum-free & Animal origin-free) is a chemically defined liquid cryopreservation medium intended for freezing and storing a variety of mammalian cell lines. This product is optimally formulated with DMSO, glucose and other nutrients, and without serum and animal origin component. The cells have high viability after resuscitation. Also it is suitable for the culture of serum-free cultured cells.

Key features

- 1.Ready-to-use serum-free freezing medium.
- 2.Suitable for long term storage of many mammalian cells.
- 3.Directly store the cells at -80°C without programmed cooling.
- 5.High cell viability after resuscitation.

Notes

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about proteins, please visit www.elabscience.com.

Focus on your research
Service for life science

1. This product is for research use only.
2. This product is sterilized by 0.22 μ m filtration.
3. It is necessary to pay attention to the aseptic operation and avoid the pollution during the culture.
4. For the freezing of stem cells (ES cells) and the cells that sensitive to DMSO, we recommend that users perform a pre-experiment before use, freeze the cells with this product for at least one week, and then perform the formal cryopreservation after confirming the performance.