

Ponceau S Solution

Cat. No: E-IR-R127

Size: 100 mL/ 200 mL/ 500 mL

Cat	Products	100 mL	200 mL	500 mL	Storage
E-IR-R127	Ponceau S Solution	100 mL	200 mL	500 mL	RT

Introduction

Ponceau S Solution can be used for protein detection on PVDF Membrane, Nitrocellulose Membrane and Cellulose Acetate Membrane, but not nylon Membrane. Ponceau S is negatively charged and can bind to positively charged amino acid residues, while Ponceau S can also bind to non-polar regions of proteins to form red bands.

This product is ready to use reagent, with low background, long-term stable performance, clear coloring and high sensitivity. More than 250 ng protein could be detected by this reagent.

Instructions

After staining, rinse with distilled water, PBS or other appropriate solution until the background is clear, and then conduct subsequent tests such as ECL, Dot-ELISA, etc.

Storage

Store at room temperature for 12 months.

Cautions

1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about products, please visit www.elabscience.com

**Focus on lab research
Service for life science**