

Ham's F-12 Medium Ingredient Formula

Ingredient(mg/L) / Cat. No.	PM150810	PM150810A	PM150811	PM150812	PM150813	PM150814
Inorganic Salt						
CaCl2 (Anhydrous)	33.3	33.3	33.3	33.3	33.3	33.3
CuSO4•5H2O	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025
FeSO4•7H2O	0.834	0.834	0.834	0.834	0.834	0.834
MgCl2 (Anhydrous)	57.22	57.22	57.22	57.22	57.22	57.22
KCl	223.6	223.6	223.6	223.6	223.6	223.6
NaHCO3	1176	1176	1176	1176	1176	1176
NaCl	7599	7599	7599	7599	7599	7100
Na2HPO4 (Anhydrous)	142	142	142	142	142	142
ZnSO4•7H2O	0.863	0.863	0.863	0.863	0.863	0.863
Amino Acid						
L-Alanyl-L-Glutamine	–	–	–	217	–	–
Glycine	7.5	7.5	7.5	7.5	7.5	7.5
L-Alanine	9	9	9	9	9	9
L-Arginine Hydrochloride	211	211	211	211	211	211
L-Asparagine Monohydrate	15.01	15.01	15.01	15.01	15.01	15.01
L-Aspartic Acid	13.3	13.3	13.3	13.3	13.3	13.3
L-Cysteine Hydrochloride Monohydrate	35.12	35.12	35.12	35.12	35.12	35.12
L-Glutamic Acid	14.7	14.7	14.7	14.7	14.7	14.7
L-Glutamine	146	146	–	–	146	146
L-Histidine Hydrochloride Monohydrate	21	21	21	21	21	21
L-Isoleucine	3.94	3.94	3.94	3.94	3.94	3.94
L-Leucine	13.1	13.1	13.1	13.1	13.1	13.1
L-Lysine Hydrochloride	36.5	36.5	36.5	36.5	36.5	36.5
L-Methionine	4.48	4.48	4.48	4.48	4.48	4.48
L-Phenylalanine	4.96	4.96	4.96	4.96	4.96	4.96
L-Proline	34.5	34.5	34.5	34.5	34.5	34.5
L-Serine	10.5	10.5	10.5	10.5	10.5	10.5
L-Threonine	11.9	11.9	11.9	11.9	11.9	11.9
L-Tryptophan	2.04	2.04	2.04	2.04	2.04	2.04
L-Tyrosine Disodium Salt Dihydrate	7.78	7.78	7.78	7.78	7.78	7.78
L-Valine	11.7	11.7	11.7	11.7	11.7	11.7
Vitamin						
Biotin	0.0073	0.0073	0.0073	0.0073	0.0073	0.0073
Choline Chloride	13.96	13.96	13.96	13.96	13.96	13.96
D-Calcium Pantothenate	0.48	0.48	0.48	0.48	0.48	0.48
Folic Acid	1.32	1.32	1.32	1.32	1.32	1.32
Niacinamide	0.037	0.037	0.037	0.037	0.037	0.037
Pyridoxine Hydrochloride	0.062	0.062	0.062	0.062	0.062	0.062
Riboflavin	0.038	0.038	0.038	0.038	0.038	0.038
Thiamine Hydrochlorid	0.34	0.34	0.34	0.34	0.34	0.34
Vitamin B12	1.36	1.36	1.36	1.36	1.36	1.36
Meso-inositol	18	18	18	18	18	18
Others						
D-Glucose	1802	1802	1802	1802	1802	1802
Thymidine	0.73	0.73	0.73	0.73	0.73	0.73
Hypoxanthine Sodium	4.77	4.77	4.77	4.77	4.77	4.77
Linoleic Acid	0.084	0.084	0.084	0.084	0.084	0.084
Lipoic Acid	0.21	0.21	0.21	0.21	0.21	0.21
Putrescine Dihydrochloride	0.161	0.161	0.161	0.161	0.161	0.161
Sodium Pyruvate	110	110	110	110	110	110
HEPES	–	–	–	–	–	5958
Phenol Red	1.2	1.2	1.2	1.2	–	1.2
Penicillin G Sodium Salt	–	100kU	–	–	–	–
Streptomycin Sulfate	–	100	–	–	–	–