

Electrophoresis Buffer (10×)

Cat. No: E-BC-R331

Size: 100 mL/ 200 mL/ 500 mL

Cat	Products	100 mL	200 mL	500 mL	Storage
E-BC-R331	Electrophoresis Buffer (10×)	100 mL	200 mL	500 mL	RT

Introduction

Tris-Glycine SDS electrophoresis buffer is used in sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE) for protein denaturation. This product is supplied with 10× concentrate for easy transportation.

Instructions

This product is a 10× concentrate, diluted with Deionized water to 1× working solution before use.
For example: take 100 mL of Electrophoresis Buffer (10×), then add 900 mL of Deionized water and mix fully.

1 × Electrophoresis Buffer Components

25 mM Tris, 192 mM Glycine, 0.1% SDS, pH8.3.

Storage

Store at room temperature for 12 months.

Cautions

1. This product may precipitate when stored at low temperature.. If precipitates, please dissolve in water bath at 37 °C.
2. Please pay attention to sealed storage to avoid contamination.
3. For your safety and health, please wear the lab coat and disposable gloves before the experiments.