



GAMMA-TYPE™



“Gamma-Type™ provides additional MHC matching information between stem cell donors and patients.”



Gamma-Type™ is a PCR SSP kit for 25 SNPs located in the Gamma genomic block of the central MHC. Gamma-Type™ is available in 25 test or 96 test kits.



Please turn over →



Gamma-Type™ will also be available with our MPS Resolver™ and Capture-HLA™ products.

Test Kit	Product Code
Gamma-Type™ 25 reaction test kit	KD PD8.0-1(25)
Gamma-Type™ 96 reaction test kit	KD PD8.0-1(96)

For Research Use Only. Not for use in diagnostic procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.

* Getz J, Sayer D, Pasquini R et al. Impact of MHC Gamma Block (Gamma Type) mismatching in the outcome of unrelated hematopoietic stem cell transplantation [abstract]. Tissue Antigens 2014; 84: 40.

Distributed by
OLERUP
www.olerup.com

Designed, developed and manufactured by Conexio.

conexio-genomics.com

MAR059 -1 Version 2.0

 **CONEXIO**
Pioneering since 2002